



PROFESSIONAL INSTRUMENT SERIES

DENTAL INSTRUMENTS CATALOGUE

2026 EDITION

Precision instruments for dental and
oral surgery workflows.



AURA SURGICARE

www.aurasurgicare.com



COMPANY PROFILE

About AURA Surgicare

A focused surgical and dental instrument catalogue developed for professional procurement and clinical supply discussions.

WHO WE ARE

AURA Surgicare supplies a broad portfolio of surgical and dental instruments for professional healthcare use. Our catalogue is structured to make instrument identification, comparison and ordering straightforward for distributors, procurement teams and clinical professionals.

OUR APPROACH

We focus on clear product references, practical instrument groupings and consistent catalogue presentation. Product selection should always be confirmed against the intended procedure, requested configuration and applicable market requirements before supply.

MANUFACTURING

AURA Surgicare supplies professional dental instruments for clinical and surgical workflows. Instrument specifications, finishes, dimensions and packaging requirements can vary by reference and order, so final technical details are confirmed during quotation and procurement.

CATALOGUE PURPOSE

This Dental Instruments Catalogue provides a structured reference for extracting forceps, elevators, diagnostic instruments, restorative instruments, orthodontic instruments, trays, hollowware and other dental instrument families.



PRODUCT WORKFLOW

From Reference to Order

A practical workflow for selecting and confirming catalogue instruments.

01

IDENTIFY

Choose the instrument family and catalogue reference that best matches the intended application.

02

CONFIRM

Confirm size, pattern, tip profile, handle style, finish and any requested configuration.

03

REVIEW

Review the requested specification and quantity before quotation and production planning.

04

INSPECT

Product and order details are checked against the confirmed specification before dispatch.

05

SUPPLY

Packaging, documentation and shipping requirements are finalized for the destination market.



CATALOGUE GUIDANCE

Using This Catalogue

Reference information is presented for identification and enquiry. Final order specifications should be confirmed before supply.

Product references

AURA reference codes are used throughout this branded edition to make enquiries and product discussions easier.

Illustrations

Instrument images and drawings are visual references. Exact appearance may vary with pattern, size, finish and production specification.

Clinical selection

Instrument selection and use remain the responsibility of appropriately trained healthcare professionals and procurement teams.

Dimensions and patterns

Where a figure, size, pattern or application is shown, use it as a catalogue reference and confirm the required configuration when ordering.

Reprocessing

Cleaning, disinfection and sterilization procedures must follow the validated instructions applicable to the supplied instrument and local requirements.



ORDERING INFORMATION

Enquiries and Product Confirmation

Use the catalogue reference together with the required quantity and configuration when requesting a quotation.

1 CATALOGUE REFERENCE

Quote the AURA product reference shown beside the instrument.

2 SIZE / PATTERN

State the required figure, size, pattern, tip or handle configuration where applicable.

3 QUANTITY

Provide the required quantity for each catalogue reference.

4 DESTINATION

Share the destination market so documentation, packaging and shipping requirements can be reviewed.


5 FINAL CONFIRMATION

Final technical and commercial details are confirmed in the quotation or order documentation.

AURA SURGICARE

www.aurasurgicare.com

TABLE OF CONTENTS




Extracting Forceps

Page 01-31




Root Elevators

Page 32-50



Filling Instruments

Page 51-55




Dental Probes

Page 56-59



Periodontia Instruments
Mallets

Page 60-67



Curettes
Scalers

Page 68-74

TABLE OF CONTENTS

Page 75-77

Modelling Instruments
Measuring Instruments




Page 78-79

Wax / Plaster Knives



Page 80-85

Diagnostic (Explorers)
Endodontic Instruments



Page 86-87

Excavators



Page 88-89

Handles For Mouth Mirrors
Mouth Mirrors



Page 90-93

Spatulas






TABLE OF CONTENTS




Needle Holders &
Saliva Ejectors

Page 94-100




Syringes & Scissors

Page 101-109



Dental Tweezers

Page 110-114




Rubber Dam Clamps

Page 115-117



Matrix Retainers &
Matrix Bands

Page 118-119




Crown Removers, Scissors,
Forceps, Amalgum Guns & Carriers

Page 120-122

TABLE OF CONTENTS

Page 123

Copper Ring Removing Pliers
& Paper Articulating Forceps




Page 124-127

Rongeurs



Page 128-129

Bone Elevators,
Files & Chisels




Page 130-151

Pliers For Orthodontics, Mouth Gags,
Tongue Depressors Lip and Cheek
Retractors and Implantology Instrumetns




Page 152-155

Impression Trays



Page 156-157

Articulators



Page 158-160

Hollowares



Page 1-31


- Extracting Forceps - English Pattern




AURA-1-5101
lower premolars



AURA-1-5102
lower premolars



AURA-1-5103
upper molars, right



AURA-1-5104
upper molars, left




AURA-1-5105
upper molars


Extracting Forceps - English Pattern



AURA-1-5106
upper wisdoms



AURA-1-5107
lower wisdoms



AURA-1-5108
lower molars



AURA-1-5109
lower molars



AURA-1-5110
Albrecht
lower molars



AURA-1-5111
lower roots



AURA-1-5112
lower roots



AURA-1-5113
lower roots



AURA-1-5114
lower roots



AURA-1-5115
lower roots, for children


Extracting Forceps - English Pattern




AURA-1-5116
upper roots




AURA-1-5117
lower roots



AURA-1-5118
upper roots



AURA-1-5119
upper roots




AURA-1-5120
upper roots



AURA-1-5121
upper roots




AURA-1-5122
upper roots



AURA-1-5123
upper roots



AURA-1-5124
upper wisdoms



AURA-1-5125
upper wisdoms


Extracting Forceps For Children - English Pattern




AURA-1-5126
Upper incisors



AURA-1-5127
Upper molars



AURA-1-5128
Upper molars, right



AURA-1-5129
Upper molars, left



AURA-1-5130
Upper roots




AURA-1-5131
lower molars



AURA-1-5132
upper incisors



AURA-1-5133
lower incisors



AURA-1-5134
upper premolars

Extracting Forceps - English Pattern




AURA-1-5135
Fractured or decayed
upper teeth



AURA-1-5136
Upper roots



AURA-1-5137
Lower roots



AURA-1-5138
Lower roots



AURA-1-5139
Lower roots



AURA-1-5140
Upper roots



AURA-1-5141
Upper roots



AURA-1-5142
Upper roots



AURA-1-5143
Upper roots



AURA-1-5144
Upper roots

Extracting Forceps - English Pattern



AURA-1-5145
Upper roots



1046
Fig 52
Upper roots and
canines




AURA-1-5147
For separating
upper molars



AURA-1-5148
For separating
upper molars



AURA-1-5149
For separating
lower molars



Extracting Forceps - English Pattern



AURA-1-5155
Upper molars, left



AURA-1-5156
Upper molars, right



AURA-1-5157
Upper third molars



AURA-1-5158
Upper third molars



AURA-1-5159
Lower roots



AURA-1-5160
Upper third molars, left




AURA-1-5161
Upper third molars, right



AURA-1-5162
Lower molars



AURA-1-5163
Lower roots



AURA-1-5164
Lower roots

Extracting Forceps - English Pattern




AURA-1-515
Lower roots




AURA-1-516
Lower canines and
premolars



AURA-1-5167
Lower third molars



AURA-1-5168
Lower third molars



AURA-1-5169
Lower third molars



AURA-1-5170
Lower third molars



AURA-1-5171
Lower third molars, left



AURA-1-5172
Lower third molars, right



AURA-1-5173
Lower molars

Extracting Forceps - English Pattern



AURA-1-5174
Lower decayed molars
or broken-down crowns




WITH SERRATED
TIPS

AURA-1-5175
Upper molars, right



WITH SERRATED
TIPS

AURA-1-5176
Upper molars, left



WITH SERRATED
TIPS

AURA-1-5177
Broken-down lower molars



WITH SERRATED
TIPS

AURA-1-5178
Upper canines



AURA-1-5179
Upper third molars



AURA-1-5180
Lower roots



AURA-1-5181
Upper roots



AURA-1-5182
Universal extracting forceps for
upper teeth, in particular molars
and premolars




AURA-1-5183
Universal extracting forceps for
lower teeth, in particular molars
and premolars



Extracting Forceps - English Pattern




AURA-1-5184
Upper canines



AURA-1-5185
Lower canines



AURA-1-5186
Lower roots



AURA-1-5187
Upper roots



AURA-1-5188
Lower roots




AURA-1-5189

Universal extracting forcep for upper incisors and canines



AURA-1-5190

Upper molars



AURA-1-5191

Lower incisors, premolar and roots

WITH SERRATED TIPS



AURA-1-5192

Lower molars

WITH SERRATED TIPS


Extracting Forceps With Anatomically Shaped Handle



AURA-1-5193
Upper incisors
and canines



AURA-1-5194
Upper incisors
and canines



AURA-1-5195
Lower incisors
and canines



AURA-1-5196
Upper premolars




AURA-1-5197
Lower premolars



AURA-1-5198
Upper molars, right



AURA-1-5199
Upper molars, left



AURA-1-5200
Lower molars



AURA-1-5201
Lower molars



AURA-1-5202
Lower roots

Extracting Forceps With Anatomically Shaped Handle



AURA-1-5203
Lower roots



AURA-1-5204
Upper roots



AURA-1-5205
Upper third molars



AURA-1-5206
Lower roots



AURA-1-5207
Lower third molars



AURA-1-5208
Lower molars



AURA-1-5209
Lower roots



AURA-1-5210
Upper molars, right




AURA-1-5211
Upper molars, left



AURA-1-5212
Lower roots



AURA-1-5213
Upper roots



AURA-1-5214
Upper roots

Extracting Forceps - American Pattern



AURA-1-5215
Lower third molars



AURA-1-5216
Separating forcep



AURA-1-5217
Upper molars



AURA-1-5218
Lower molars



AURA-1-5219
Lower molars



AURA-1-5220
Upper molars, left



AURA-1-5221
Upper molars, right



AURA-1-5222
Lower molars



AURA-1-5223
Upper molars, right



AURA-1-5224
Upper molars, left

Extracting Forceps - American Pattern



AURA-1-5225
Upper incisors and roots



AURA-1-5226
Upper and lower roots



AURA-1-5227
Upper molars, right



AURA-1-5228
Upper molars, left



AURA-1-5229
Deciduous teeth



AURA-1-5230
Upper incisors, premolars, roots



AURA-1-5231
Upper teeth (for children)



AURA-1-5232
Upper incisors, premolars, roots



AURA-1-5233
Lower teeth (for children)




AURA-1-5234
Lower third molars

Extracting Forceps - American Pattern




AURA-1-5235
upper incisors and bicuspids



AURA-1-5236
upper roots, anterior



AURA-1-5237
lower wisdoms



AURA-1-5238
lower molars, splitting forceps



AURA-1-5239
upper molars



AURA-1-5240
Lower premolars

WITH SERRATED
TIPS



AURA-1-5241
Lower molars



AURA-1-5242
Upper incisors
and roots



AURA-1-5243
Lower roots

WITH SERRATED
TIPS



AURA-1-5244
Lower incisors
and canines

WITH SERRATED
TIPS

Extracting Forceps For Children - Klein Pattern



AURA-1-5245
Upper incisors



AURA-1-5246
Upper premolars



AURA-1-5247
Upper molars



AURA-1-5248
Upper roots



AURA-1-5249
Lower incisors



AURA-1-5250
Lower molars



AURA-1-5251
Lower roots



AURA-1-5252
Lower roots



AURA-1-5253
Upper roots



AURA-1-5255
Roots fragments




AURA-1-5254
Roots fragments

Page 32–50

- Root Elevators

lower jaw, left side

lower jaw, right side



Root elevator, spiral-shaped, for extremely facilitated extraction of the molars in the lower jaw. The slender, conical points allow deep penetration between tooth and gum. Owing to the hook-shaped working parts, slipping of the elevator is avoided.



AURA-2-5101



AURA-2-5102

Root Elevators



Bein



AURA-2-5103



AURA-2-5104



AURA-2-5105



AURA-2-5106



AURA-2-5107



AURA-2-5108



AURA-2-5109



Seldin



AURA-2-5110



AURA-2-5111



AURA-2-5112



AURA-2-5113



AURA-2-5114



AURA-2-5115

Root Elevators



AURA-2-5116



AURA-2-5117



AURA-2-5118



AURA-2-5119



AURA-2-5120



AURA-2-5121



AURA-2-5122



Apical



AURA-2-5123



AURA-2-5124



AURA-2-5125



AURA-2-5126

Root Elevators



Cryer



AURA-2-5127



AURA-2-5128



Coupland



AURA-2-5129



AURA-2-5130



AURA-2-5131



AURA-2-5132



AURA-2-5133



AURA-2-5134



AURA-2-5135



AURA-2-5136



AURA-2-5137

Root Elevators



AURA-2-5138



AURA-2-5139



AURA-2-5144



AURA-2-5145



Kopp



AURA-2-5140



AURA-2-5141



Kopp



AURA-2-5146



AURA-2-5147



AURA-2-5142




AURA-2-5143



AURA-2-5148



AURA-2-5150



AURA-2-5151




Hylin



AURA-2-5149



Bein



AURA-2-5152



AURA-2-5153

Root Elevators



mes
AURA-2-5154



mes
AURA-2-5155



mes
AURA-2-5156



London
AURA-2-5157




London
AURA-2-5158



London
AURA-2-5159



AURA-2-5160



Pott



AURA-2-5162

AURA-2-5163



AURA-2-5164



AURA-2-5165

Root Elevators



AURA-2-5166



AURA-2-5167




AURA-2-5168



AURA-2-5169



AURA-2-5170



AURA-2-5171



AURA-2-5172



AURA-2-5173



AURA-2-5174



AURA-2-5175



AURA-2-5176

Root Elevators Root Tip Elevators



AURA-2-5177



AURA-2-5178



AURA-2-5179



AURA-2-5180



AURA-2-5181



AURA-2-5182



AURA-2-5183




AURA-2-5184



AURA-2-5185



AURA-2-5186



Winter



AURA-2-5187



AURA-2-5188



AURA-2-5189



AURA-2-5190



AURA-2-5191



AURA-2-5192

Root Elevators



AURA-2-5193



AURA-2-5194



AURA-2-5195



AURA-2-5196



Barry



AURA-2-5197



AURA-2-5198



AURA-2-5199



AURA-2-5200



AURA-2-5201



AURA-2-5202



AURA-2-5203



AURA-2-5204



AURA-2-5205



AURA-2-5206



AURA-2-5207

Root Elevators



AURA-2-5208



AURA-2-5209



AURA-2-5210



AURA-2-5211



AURA-2-5212



AURA-2-5213



AURA-2-5214



AURA-2-5215



AURA-2-5216



AURA-2-5217



AURA-2-5218



AURA-2-5219



AURA-2-5220



AURA-2-5221

Root Elevators

Root Elevators according to Dr. Sader,
with interchangeable working tips



AURA-2-5222




AURA-2-5223




AURA-2-5224

Page 51–55

- Filling Instruments



AURA-2-5226 AURA-2-5227 AURA-2-5228 AURA-2-5229




AURA-2-5234 AURA-2-5235 AURA-2-5236



AURA-2-5225


AURA-2-5230 AURA-2-5231 AURA-2-5232




AURA-2-5233


AURA-2-5237 AURA-2-5238 AURA-2-5239

Filling Instruments





Filling Instruments



AURA-2-5268
Fig 1



AURA-2-5269
Fig 2



AURA-2-5270
Fig 3



AURA-2-5271



AURA-2-5272




AURA-2-5274



AURA-2-5273




AURA-2-5275



AURA-2-5276



AURA-2-5277



AURA-2-5278
Fig 160



AURA-2-5279



AURA-2-5280

Fig 1



AURA-2-5281

Fig 2



AURA-2-5282

Fig 3



AURA-2-5283

Fig 4



AURA-2-5284

Fig 5



AURA-2-5285

Fig 6



AURA-2-5286



AURA-2-5287



AURA-2-5288



AURA-2-5289



AURA-2-5290



AURA-2-5291



AURA-2-5292



AURA-2-5293




AURA-2-5294

Page 56–59

- Dental Probes

Dental Probes

With Point Set In



AURA-2-5295



AURA-2-5296



AURA-2-5297



AURA-2-5298



AURA-2-5299



AURA-2-5300



AURA-2-5301



AURA-2-5302



AURA-2-5303



AURA-2-5304



AURA-2-5305




AURA-2-5306



AURA-2-5307

Polishing Instruments



Periodontal Pocket Probes



- New marking process.
- Laser system marking.
- No more paint peeling off.
- Permanent marking.
- Antiglare, low depth marking.
- Maximum tip resistance against breaking.



AURA-2-5322



AURA-2-5323



AURA-2-5324



AURA-2-5325

O.M.S



AURA-2-5326



AURA-2-5327



AURA-2-5328



AURA-2-5329
Fox



AURA-2-5330
Michigan



AURA-2-5331

Page 60–67

- Periodontia Instruments, Mallets



Files

Chisels



AURA-3-5101
Fig. 1S/2S



AURA-3-5102
Fig. 3S/4S



AURA-3-5103
Fig. 9/10



AURA-3-5104
Fig. 11/12



AURA-3-5105
Fig. 3



AURA-3-5106
Fig. 4



AURA-3-5107
Back - Action
Fig. 36/37

Periodontia Instruments

Files



AURA-3-5108
Fig. 1S/2S



AURA-3-5109
Fig. 3S/4S



AURA-3-5110
Fig. 9/10



AURA-3-5111
Fig. 11/12



AURA-3-5112
Fig. 3



AURA-3-5113
Fig. 4



AURA-3-5114
Back - Action
Fig. 36/37

Gingivectomy knives



AURA-3-5115
Fig. 7



AURA-3-5116
Fig. 8



AURA-3-5117
Fig. 9



AURA-3-5118
Fig. 11



AURA-3-5119
Fig. 3/4



AURA-3-5120
Fig. 5/6

Periodontia Instruments

Curettes and Scalers



AURA-3-5121
Fig. HU 1/2



AURA-3-5122
Fig. 2R/2L



AURA-3-5123
Fig. 4R/4L



AURA-3-5124
Fig. 13/14



AURA-3-5125



AURA-3-5126



AURA-3-5127
el

Gingivectomy knives



AURA-3-5128
ox
Fig. 7



AURA-3-5129
ox
Fig. 8



AURA-3-5130
ox
Fig. 9



AURA-3-5131
ox
Fig. 11



AURA-3-5132
Fig. 3/4



AURA-3-5133
Fig. 5/6

Periodontia Instruments

Curettes and Scalers



AURA-3-5134
Fig. HU 1/2



AURA-3-5135
Fig. 2R/2L



AURA-3-5136
Fig. 4R/4L



AURA-3-5137
Fig. 13/14



AURA-3-5138



AURA-3-5139



AURA-3-5140
el


Bone Chisels and Gouges




AURA-3-5141	AURA-3-5142	AURA-3-5143	AURA-3-5144	AURA-3-5145	AURA-3-5146	AURA-3-5147	AURA-3-5148	AURA-3-5149	AURA-3-5150
Fig. 1	Fig. 2	Fig. 1	Fig. 2	Fig. 5	Fig. 5R	Fig. 6	Fig. 1W	Fig. 2W	Fig. 3

Periodontia Instruments and Mallets

Bone Chisels and Gouges




Martelli chirurgici Mallet



Page 68–74

▪ Currettes Scalers





AURA-3-5170
Fig. 1/2



AURA-3-5171
Fig. 3/4



AURA-3-5172
Fig. 5/6



AURA-3-5173
Fig. 7/8



AURA-3-5174
Fig. 9/10



AURA-3-5175
Fig. 11/12



AURA-3-5176
Fig. 13/14




AURA-3-5177
Mesial Distal
Fig. 11/14



AURA-3-5178
Mesial Distal
Fig. 12/13

Gracey Curettes AF Type



AURA-3-5179
AF/type
Fig. 1/2



AURA-3-5180
AF/type
Fig. 3/4



AURA-3-5181
AF/type
Fig. 5/6



AURA-3-5182
AF/type
Fig. 7/8




AURA-3-5183
AF/type
Fig. 11/12




AURA-3-5184
AF/type
Fig. 13/14



AURA-3-5185
AF/type
Fig. 15/16




AURA-3-5187	AURA-3-5188	AURA-3-5189	AURA-3-5190	AURA-3-5191	AURA-3-5192	AURA-3-5193
AF/type	AF/type	AF/type	AF/type	AF/type	AF/type	AF/type
Fig. 1/2	Fig. 3/4	Fig. 5/6	Fig. 7/8	Fig. 11/12	Fig. 13/14	Fig. 15/16



AURA-3-5186							
AF/type							
	AURA-3-5194	AURA-3-5195	AURA-3-5196	AURA-3-5197	AURA-3-5198	AURA-3-5199	AURA-3-5200
	MF/type	MF/type	MF/type	MF/type	MF/type	MF/type	MF/type
	Fig. 1/2	Fig. 3/4	Fig. 5/6	Fig. 7/8	Fig. 11/12	Fig. 13/14	Fig. 15/16

Gracey Curettes AF Type



AURA-3-5201
MF/type
Fig. 1/2



AURA-3-5202
MF/type
Fig. 3/4



AURA-3-5203
MF/type
Fig. 5/6



AURA-3-5204
MF/type
Fig. 7/8



AURA-3-5205
MF/type
Fig. 11/12



AURA-3-5206
MF/type
Fig. 13/14



AURA-3-5207
MF/type
Fig. 15/16



AURA-3-5208_{ax}
Fig. 1



AURA-3-5209_{ax}
Fig. 2




AURA-3-5210_{ax}
Fig. 3



AURA-3-5211_{ax}
Fig. 4




AURA-3-5212_{ax}
Fig. 21



AURA-3-5213_{ax}
Fig. 5




AURA-3-5214
Fig. 13S/14S



AURA-3-5215_{od}
Fig. 7/8

Curettes and Scalers





AURA-3-5228
Fig. 1



AURA-3-5229
Fig. 2



AURA-3-5230
Fig. 3



AURA-3-5231
Fig. 85



AURA-3-5232
Fig. 86




AURA-3-5233
Fig. 87


Page 75–77

- Modelling Instruments & Measuring Instruments





AURA-4-5101 AURA-4-5102 AURA-4-5103 AURA-4-5104 AURA-4-5105 AURA-4-5106 AURA-4-5107 AURA-4-5108 AURA-4-5109 AURA-4-5110 AURA-4-5111 AURA-4-5112



AURA-4-5113 AURA-4-5114 AURA-4-5115 AURA-4-5116 AURA-4-5117 Fig, 1 AURA-4-5118 Fig, 2 AURA-4-5119 Fig, 3 AURA-4-5120 Fig, 4 AURA-4-5121 Fig, 5

Measuring Instruments



AURA-4-5122




AURA-4-5123
130 mm



AURA-4-5124
12 cm



AURA-4-5125




AURA-4-5126



AURA-4-5127
Scleral Marker



AURA-4-5128



AURA-4-5129



AURA-4-5130
7 1/8"



AURA-4-5131
7 1/8"



AURA-4-5132
6 3/4"
Gritman



AURA-4-5133
6 3/4"
Medicon



AURA-4-5134
8 3/4"



AURA-4-5135



AURA-4-5136



AURA-4-5137

Page 78–79

▪ Wax Knives/plaster



Wax Knives



AURA-4-5138_k
Fig. 1
12,5cm



AURA-4-5139_k
Fig. 2
17cm



AURA-4-5140
Fig. 1
12,5cm



AURA-4-5141
Fig. 2
17cm



AURA-4-5142_{ino}
Fig. 1
12,5cm



AURA-4-5143_{ino}
Fig. 2
17cm



AURA-4-5144
17 cm



AURA-4-5145
Fig, 2
17 cm



AURA-4-5146
16 cm



AURA-4-5147
16 cm


Page 80–85

▪ Diagnostic, Explorer, Endodontic Instruments



Diagnostic

Explorers



AURA-4-5148

AURA-4-5149

AURA-4-5150

AURA-4-5151


AURA-4-5152

AURA-4-5153

AURA-4-5154

AURA-4-5155

Rein



AURA-4-5156
Spreader

AURA-4-5157
Spreader

AURA-4-5158
Spreader

AURA-4-5159
Spreader

AURA-4-5110
Spreader


AURA-4-5111
Spreader

AURA-4-5112
Spreader

Endodontic Instruments

Double-Ended Root Canal Explorers Plugger/Spreaders

Ø 0.4 Ø 0.4 Ø 0.6



Spreaders


Ø 0.6 Ø 0.5



Periodontal Pocket Probes



Root-canal pluggers



Endodontic set




AURA-4-5132
late


Endodontic Instruments

Root-canal pluggers

Optimum Flexibility
Very Fine Points



Root-canal explorers



AURA-4-5143

AURA-4-5144

AURA-4-5145

Root-canal excavators



AURA-4-5146

AURA-4-5147

AURA-4-5148

Page 86–87

- Excavators



Excavators

1 mm



AURA-4-5149

1.2 mm



AURA-4-5150

1 mm



AURA-4-5151

1.2 mm



AURA-4-5152

1.6 mm



AURA-4-5153

2 mm



AURA-4-5154

1 mm



AURA-4-5155
Fig.243

1,4 mm



AURA-4-5156
Fig.G2

0,8 mm



AURA-4-5157

0,8 mm



AURA-4-5158

1 mm



AURA-4-5159

1 mm



AURA-4-5160

1,2 mm



AURA-4-5161

1,2 mm



AURA-4-5162

1,6 mm



AURA-4-5163

1,6 mm



AURA-4-5164

2 mm



AURA-4-5165

2 mm




AURA-4-5166

Page 88–89

- Handles For Mouth Mirrors





Handles for Mouth Mirrors and Mouth Mirrors



AURA-5-5112



AURA-5-5113



AURA-5-5114




AURA-5-5115




AURA-5-5116



AURA-5-5117^{ter}
mouth mirror



AURA-5-5118^{ter}
mouth mirror
+ handle
code 2101



AURA-5-5123^{et}



AURA-5-5124^{et}




AURA-5-5125^{et}

Page 90–93

- Spatulas





AURA-5-5126




AURA-5-5127



AURA-5-5128



AURA-5-5129



AURA-5-5130



AURA-5-5131





AURA-5-5132




AURA-5-5133

Cement Spatulas





Spatulas for Plaster and Alginate



Page 94–100

- Needle Holders





AURA-5-5150
16 cm



AURA-5-5151
16 cm



AURA-5-5152
18 cm



AURA-5-5153
18 cm



AURA-5-5154
18 cm



AURA-5-5155
15 cm



AURA-5-5156
12 cm



AURA-5-5157
14 cm

Needle Holders



AURA-5-5158
15 cm



AURA-5-5159
TC
15 cm



AURA-5-5160
17.5 cm



AURA-5-5161
r
14 cm



AURA-5-5162
r
14.5 cm



AURA-5-5163
16 cm



AURA-5-5164
16 cm



AURA-5-5165 Hegar-
r 14 cm



AURA-5-5166
18 cm



AURA-5-5167
18 cm



AURA-5-5168
17 cm



AURA-5-5169
17 cm



AURA-5-5170
19.5 cm
athieu



AURA-5-5171
12 cm



AURA-5-5172
16 cm



AURA-5-5173
12 cm



AURA-5-5174
16 cm

Needle Holders



AURA-5-5175
14 cm



AURA-5-5176
14 cm



AURA-5-5177
14/60 cm



AURA-5-5178
18 cm



AURA-5-5179
18 cm



AURA-5-5180
14 cm



AURA-5-5181
14 cm



AURA-5-5182
14/60 cm



AURA-5-5183
18 cm



AURA-5-5184
18 cm



AURA-5-5185
14 cm



AURA-5-5186
14 cm



AURA-5-5187
14 cm



AURA-5-5188
14 cm



AURA-5-5189
17 cm



AURA-5-5190
17 cm



AURA-5-5191
14 cm



AURA-5-5192
18 cm
quer

Needle Holders & Stainless Steel Saliva Ejectors



The tips have tungsten carbide inserts for better handling.



AURA-5-5193
17 cm



AURA-5-5194
17 cm



AURA-5-5195
14 cm



AURA-5-5196
13 cm



AURA-5-5197 Barraquer
13 cm



AURA-5-5198
13 cm



AURA-5-5199
13 cm



AURA-5-5200



AURA-5-5201




AURA-5-5202




AURA-5-5206
Fraizer
15.02mm,cm

0.4 AURA-5-5207




Derf

0.4 AURA-5-5208




Halsey

0.5 AURA-5-5209




Baumgartner

0.4 AURA-5-5210




Crilewood

0.4 AURA-5-5211



Swedish-Micro
Vascular

0.5 AURA-5-5212



Mayo-Hegar


0.5 AURA-5-5213

0.4 AURA-5-5214



Mathieu
Kocher

0.4 AURA-5-5215




Mathieu

0.4 AURA-5-5216



0.4 AURA-5-5217




Mathieu


0.4 AURA-5-5218



0.4 AURA-5-5219



0.4 AURA-5-5220



Olsen-Hegar

0.5 AURA-5-5221



Toennis

Page 101–109

- Syringes and Scissors



Water Syringes



AURA-5-5222
Water Syringe complete with
1 cannula (Lure-Lock)



AURA-5-5223
Water Syringe complete with
1 cannula (Lure-Lock)



AURA-5-5224
Metal water syringe, self-filling,
with 1 cannula (Lure-Lock)



AURA-5-5225
Sparrel Barrel with
plunger suitable for
syringes

AURA-5-5226
Cannula,
only, with
Luer-Lock connector
suitable for syringes



AURA-5-5227
Syringe



AURA-5-5228
Band Adaptor



AURA-5-5229
Syringe



AURA-5-5230
Syringe



AURA-5-5231
Ruler



AURA-5-5232
Ruler
With Adjustable Slide



AURA-5-5233
Syringe

Dental Syringes




AURA-5-5234



AURA-5-5235



AURA-5-5236



AURA-5-5237 1,8 ml
AURA-5-5238 2,2 ml
Cartridge Syringe complete
With 2 tips and without needles



AURA-5-5239
Tip Only, short



AURA-5-5240
Tip Only, long



AURA-5-5241



AURA-5-5242



AURA-5-5243



AURA-5-5244
tip



AURA-5-5245
tip




AURA-5-5246
screw tip , cork



AURA-5-5247
screw tip , cork



AURA-5-5248









Operating Scissors



Suture Scissors | Tijeras para Hilos | Fadenscheren







Scissors

Scissors, One Blade Serrated

- One blade is with fine serrations to avoid tissue or ligature slippage.
- Manufactured from high quality stainless steel to provide longer-lasting cutting edge
- Single black ring for easy identification.



Iris
11.5 cm




Goldman-Fox
13 cm

AURA-5-5305




AURA-5-5306
13 cm



AURA-5-5307



AURA-5-5308



Kelly
16 cm




AURA-5-5309
16 cm



AURA-5-5310
16 cm

Super-Cut

- Cutting Blade horned to "razor" sharpness with opposing blade finished with fine serrations to avoid tissue slippage.
- Extremely smooth cutting with less trauma.
- Manufactured from high quality stainless steel to provide longer-lasting cutting edge.
- Single gold plated ring for easy identification.




AURA-5-5311



AURA-5-5312



Iris
11.5 cm



AURA-5-5313
11.5 cm









AURA-5-5314



AURA-5-5315




Goldman-Fox
13 cm



Micro Scissors



Micro Scissor
Stainless Steel




AURA-5-5327

AURA-5-5328



Micro Scissor
Stainless Steel




AURA-5-5329

AURA-5-5330




Micro Scissor
Stainless Steel
for ligature, one blade serrated




AURA-5-5331


Micro Scissors with Titanium Stainless Steel Tips




AURA-5-5332




AURA-5-5333



AURA-5-5334



AURA-5-5335



AURA-5-5336

Page 110–114

- Dental Tweezers





AURA-5-5337
15 cm



smooth



AURA-5-5338
15 cm



serrated



AURA-5-5339
15 cm



smooth



AURA-5-5340
15 cm



serrated



AURA-5-5341
15 cm



smooth



AURA-5-5342
15 cm



serrated



AURA-5-5343
15 cm



smooth



AURA-5-5344
15 cm



serrated

Dental Tweezers



AURA-5-5345
16 cm



smooth



AURA-5-5346
16 cm



serrated



AURA-5-5347
16 cm



smooth



AURA-5-5348
16 cm



serrated



AURA-5-5349
16 cm



smooth



AURA-5-5350
16 cm



serrated



AURA-5-5351
16 cm




smooth



AURA-5-5352
16 cm



serrated




AURA-5-5353
14.5 cm, 5³/₄"

AURA-5-5354
16 cm, 6¹/₄"

AURA-5-5355
16 cm, 6¹/₄"

AURA-5-5356
15 cm, 6"

AURA-5-5357
1¹/₈"



AURA-5-5358
1¹/₈"

AURA-5-5359
15 cm, 6" lege

AURA-5-5360
15 cm, 6" lege

AURA-5-5361
15 cm, 6" lege

AURA-5-5362
15 cm, 6" lege

Dental Tweezers



AURA-5-5363
15 cm, 6"



AURA-5-5364
15 cm, 6"



AURA-5-5365
16 cm, 6 1/4"



AURA-5-5366
16 cm, 6 1/4"



AURA-5-5367
15.5 cm, 6 1/8"



AURA-5-5368
15.5 cm, 6 1/8"



AURA-5-5369



AURA-5-5370
1/2



AURA-5-5371
1/8



AURA-5-5372
1/8



AURA-5-5373
15 cm, 6"




AURA-5-5374
15 cm, 6"



AURA-5-5375
16 cm, 6 1/4"



AURA-5-5376
18 cm, 7 1/8"



AURA-5-5377
16 cm, 6 1/4"



AURA-5-5378
16 cm, 6 1/4"



AURA-5-5379
1/4"



AURA-5-5380
1/4"




AURA-5-5381
1/4"

Page 115–117


- Rubber Dam Clamps



Rubber Dam Clamp Forceps



CLAMPS (PLIERS)





AURA-5-5392
Labial Clamp



AURA-5-5393 Winged al Clamp



AURA-5-5394
Cervical Clamp



AURA-5-5395
Labial Clamp

Premolars



AURA-5-5396



AURA-5-5397



AURA-5-5398

For Small Teeth, Elongated or With Narrow Collars



AURA-5-5399



AURA-5-5400



AURA-5-5401

For General Use With Adjustable Jaws According To Gingival Line Shapes



AURA-5-5402



AURA-5-5403
ged



AURA-5-5404

For Premolars, With Elongated Collar



AURA-5-5405



AURA-5-5406
ged



AURA-5-5407

For Premolars, With Elongated Collar



AURA-5-5408



AURA-5-5409

Rubber Dam Clamps



AURA-5-5410

AURA-5-5411

AURA-5-5412

AURA-5-5413

For Small Lower Molars

For Small Lower Molars With Round clamps



AURA-5-5414

AURA-5-5415

AURA-5-5416

AURA-5-5417

For Square Lower Molars

For Round Upper Molars



AURA-5-5418

AURA-5-5419

AURA-5-5420

AURA-5-5421

For Molar Roots

For Partially Broken Molars



AURA-5-5422

AURA-5-5423
ed

AURA-5-5424

AURA-5-5425

AURA-5-5426

For Round Molars, Partially Broken, Undersized or Inclined



AURA-5-5427

AURA-5-5428

AURA-5-5429

AURA-5-5430


Page 118–122

- Matrix Retainer
- Matrix Bands





AURA-6-5101
63mm




AURA-6-5102
universal



AURA-6-5103
unior



AURA-6-5104
enior



AURA-6-5105
5 mm



AURA-6-5106
6 mm



AURA-6-5107
7 mm



AURA-6-5108



AURA-6-5109



AURA-6-5111



AURA-6-5112



AURA-6-5113



AURA-6-5114

1750-5 mm x 100 cm
1750-6 mm x 100 cm
1750-7 mm x 100 cm



AURA-6-5110

For Molars and Premolars



AURA-6-5115



AURA-6-5116



AURA-6-5117

Matrix Bands



AURA-6-5118

Fig 7



AURA-6-5119

Fig 8



AURA-6-5120

Fig 9



AURA-6-5121

Fig 10



AURA-6-5122



AURA-6-5123



AURA-6-5124



AURA-6-5125



AURA-6-5126



AURA-6-5127



AURA-6-5128



AURA-6-5129



AURA-6-5130



AURA-6-5134



AURA-6-5131



AURA-6-5132





AURA-6-5133

Page 120–122

- Crown Removers
- Scissors
- Forceps
- Amalgum Guns
- Carries





Dental Crown Removers, Scissors, Forceps, Amalgam Guns & Carriers



Universal "HM"
12 cm, 4³/₄"



AURA-6-5144



AURA-6-5145



AURA-6-5146
10,5 cm, 4¹/₈"



AURA-6-5147



AURA-6-5148



AURA-6-5149




AURA-6-5150



AURA-6-5151



AURA-6-5152




AURA-6-5153
Crown Remover
With 2 Points




AURA-6-5154



AURA-6-5155



AURA-6-5156 ge 20 cm	AURA-6-5157 bo 20 cm	AURA-6-5158 o 20 cm	AURA-6-5159 bo 20 cm	Ø 2 mm AURA-6-5160 15.5 cm	Ø 2.5 mm AURA-6-5161 15.5 cm	Ø 3 mm AURA-6-5162 15.5 cm
----------------------------	----------------------------	---------------------------	----------------------------	----------------------------------	------------------------------------	----------------------------------



Ø 2 mm AURA-6-5163 15 cm	Ø 2 mm AURA-6-5164	AURA-6-5165 90°
--------------------------------	-----------------------	--------------------



AURA-6-5166
ell
Non Skid Base
Ø 3.5 mm

Page 123

- Copper Ring Removing Plierspaper Articulator Forceps

Copper Ring Removing Pliers, Paper Articulator Forceps



AURA-6-5167
11,5 cm



AURA-6-5168
14,5 cm

Copper ring removing pliers



AURA-6-5169
16 cm



AURA-6-5170
16 cm



AURA-6-5171
15 cm



AURA-6-5172
15 cm



AURA-6-5173
14 cm




AURA-6-5174
12 cm osquito




AURA-6-5175
10 cm uito

Page 124-127


- Rongeurs




AURA-6-5176
14 cm, 5 1/2"




AURA-6-5177
1 1/2"




AURA-6-5178
14 cm, 5 1/2"



AURA-6-5179
14 cm, 5 1/2"



AURA-6-5180
17 cm, 6 3/4"



AURA-6-5181
17 cm, 6 3/4"

Bone Rongeurs



AURA-6-5182
15 cm



AURA-6-5183
14 cm



AURA-6-5184
12 cm



AURA-6-5185
10 cm



AURA-6-5186
15 cm



AURA-6-5187
14,5 cm



AURA-6-5188
14 cm



AURA-6-5189
15 cm




AURA-6-5190
15 cm




AURA-6-5191
16 cm, 6 1/4"



AURA-6-5192
16 cm, 6 1/4"



AURA-6-5193
16 cm, 6 1/4"



AURA-6-5194_Z
14 cm, 5 1/2"



AURA-6-5195
15 cm



2 mm



AURA-6-5196
18 cm, 7 1/8"




3 mm



AURA-6-5197_{en}
18 cm, 7 1/8"


Rongeurs




AURA-6-5198
No. 5 S, 14 cm, 5 1/2"




AURA-6-5199
No. 5 S, 17 cm, 6 3/4"




AURA-6-5200
16 cm, 6"



AURA-6-5201
16 cm, 6 1/4"



AURA-6-5202
16 cm, 6 1/4"




AURA-6-5203
15 cm, 6"

Page 128–129


- Bone Elevators, Bone Files and Chisels





AURA-6-5204
16 cm, 6 1/4"

AURA-6-5205
17 cm, 6 3/4"





AURA-6-5206
19 cm, 7 1/2"

AURA-6-5207
17 cm, 6 3/4"

AURA-6-5208
19 cm, 7 1/2"

AURA-6-5209
19,5 cm, 7 3/4"

Bone Files and Chisels



Page 130–151

- Pliers For Orthodontics
- Mouth Gags
- Tongue Depressors
- Lip and Cheek Retractors
- Implantology Instruments



Pinze per formare loops
Loops forming pliers



AURA-7-5101

Adams pliers for bending
wires up to 0,7mm. (.28\"



AURA-7-5102

Adams pliers for bending
wires up to 0,7mm. (.28\"



AURA-7-5103

Tweed loop forming pliers.
For forming loop in wires up
to 0,56 mm. (.22\"



AURA-7-5104

Aderer three jaws bending
pliers. For wires up
to 0,7 mm. (.28\"



AURA-7-5105

Aderer three jaws bending
pliers. For wires up to 0,7 mm.
(.28\"

Pliers for Othodontics

Ligature Cutter Pliers



AURA-7-5106

Micro-miniature light wire cutter. It cuts soft ligature wires. For use in caso of limited access. Caution: maximum cutting capacity: soft ligature wire up to 0,30 mm. (.012").



AURA-7-5107

Wire cutter. For ligature wires and pins. It is strong enough to bend archwire ends after distal cutting. Maximum cutting capacity up to 0,56 x 0,7 mm. (.022" x .028").



AURA-7-5108

Wire cutter. For ligature wires and pins. It is strong enough to bend archwire ends after distal cutting. Maximum cutting capacity: up to 0,9 mm. (.036").



AURA-7-5109

Wire cutter. For ligature wires and pins. It is strong enough to bend archwire ends after distal cutting. Maximum cutting capacity: up to 1,5 mm. (.044").

Pinze per formare laboratorio
Laboratory pliers



AURA-7-5110

Flat nose pliers with smooth beaks for bending wires up to 0,7 mm. (.028").



AURA-7-5111

Bending wires up to 0,7 mm. (.028").



AURA-7-5112


Bending wires up to 0,7 mm. (.028").



AURA-7-5113

Bending wires up to 0,7 mm. (.028").

Pliers for Othodontics



AURA-7-5114

Nance loop closing pliers, for bending of stop loops and spe ms' clasps. Wires uo to 0,7 mm. (.028").



AURA-7-5115

Universal pliers for bending wires up to 0,9 mm. (.035") or cutting up to 0,7 mm. (.028").



AURA-7-5116

Stainless steel pliers tray standard dimensions



AURA-7-5117

Universal pin roughening, cutting and bending pliers. Forwires uo to 0,9 mm. (.035"). Will cut wires up to 0,7 mm. (.028").

Utility and specialty pliers



AURA-7-5118

Reynolds contouring pliers.



AURA-7-5119

Lingual arch pliers used to form double and triple back bends in 0,7 mm. (.028") or 0,9 mm. (.035") wire for insertion in lingual sheaths.



AURA-7-5120

Weingart pliers with serrated points.



AURA-7-5121

De La Rosa arch forming pliers. Designed to form and contour arch wires on round or rectangular wire. The beak has grooves of 0,4 mm (.016"), 0,46 mm (.018") and 0,56 mm (.022").



AURA-7-5122

Arch forming pliers with a non grooved contouring surface the forming of precise arches.

Pliers for Othodontics



AURA-7-5123

Waldsachs universal pliers for bending wires up to 1,0 mm. (.039") or cutting up to 0,7 mm. (.028").



AURA-7-5124

Round and Concave pliers for loop forming and bending of thin round and square wires up to 0,5 mm. (.020").

Loop formare pliers



AURA-7-5125

Tweed loop Forming Pliers. For round or rectangular wires up to 0,56 x 0,7 mm. (.022" x .028"). Replaceable cylindrical beak with three sizes: 1,15 mm. (.045), 1,50 mm. (.059"), 1,90 mm. (.075"). Concave beak with slight parallel serrations, to keep the wire in proper position.



AURA-7-5126

Tweed Loop Forming Pliers. For making perfect Omega Loops with round or rectangular wires up to 0,56 x 0,7 mm. (.022" x .028"). The concave beak is exactly the same as pliers 5809-1



AURA-7-5127

Single diameter loop forming pliers. Wire up to 0,56 x 0,7 mm. (.022" x 0.28"). Replaceable concave beak has only one diameter section of 1,50 mm. (.059").

Distal end cutters



AURA-7-5128

Distal end cutter -safety hold- a precision cutter designed for distal arch cutting in the patient's mouth. It cuts hard wires end holds the loose distal end preventing any damaging of oral tissue. Maximum cutting capacity rectangular wires up to 0,56 x 0,7 mm. (.022" x .020").



AURA-7-5129
Micro- distal hold-. Similar to 3227 safety



AURA-7-5130
Direct bonding er (wide). The removal grip firmly for easy incisal-lingival (vertical) use or, if needed, for mesial-distal (horizontal) use with the same ease.



AURA-7-5131

Direct bond attachment remover (angled). For removing all bonded attachments. May be used mesial-distal-horizontal) or incisal-lingival (vertical) in anterior areas. Angled perfectly for easy access to difficult posterior areas.

Pliers for Othodontics

Bird beack pliers



AURA-7-5132
Wire bending pliers. For
wires up to 0,9 mm. (.036").




AURA-7-5133
Bird Beak Angle pliers for bending
round, square or rectangular wire
up to 0,7 mm. (.028").



AURA-7-5134
Double rounded jaw wire
bending pliers. For wires up
to 0,7 mm.(.028").



AURA-7-5135
Bird beak pliers for bending round, square or
rectangular wire up to 0,7 mm. (.028"). With
cutter.



AURA-7-5136
Double rounded jaw wire
bending pliers. For wires
up to 0,7 mm. (.028").



AURA-7-5137
125 mm



AURA-7-5138
125 mm



AURA-7-5139
140 mm



AURA-7-5140
145 mm



AURA-7-5141
125 mm



AURA-7-5142
125 mm



AURA-7-5143
130 mm



AURA-7-5144
130 mm



AURA-7-5145
140 mm



AURA-7-5146
130 mm




AURA-7-5147
155 mm




AURA-7-5148
150 mm



AURA-7-5149
155 mm



AURA-7-5150
125 mm



AURA-7-5151
130 mm

Pliers for Othodontics and Prosthetics



AURA-7-5152
140 mm



AURA-7-5153
130 mm



AURA-7-5154
130 mm



AURA-7-5155
125 mm



AURA-7-5156
125 mm



AURA-7-5157
140 mm



AURA-7-5158
130 mm



AURA-7-5159
130 mm



AURA-7-5160
130 mm



AURA-7-5161
125 mm



AURA-7-5162



AURA-7-5163



AURA-7-5164



AURA-7-5165



AURA-7-5166




AURA-7-5167




AURA-7-5168




AURA-7-5169




AURA-7-5170




AURA-7-5171



AURA-7-5172



AURA-7-5173



AURA-7-5174



AURA-7-5175



AURA-7-5176



AURA-7-5177




AURA-7-5178



AURA-7-5179



AURA-7-5180



AURA-7-5181



AURA-7-5182



AURA-7-5183



AURA-7-5184



AURA-7-5185



AURA-7-5186

Mouth Gags



AURA-7-5187 10 cm
AURA-7-5188 11 cm
AURA-7-5189 14 cm



AURA-7-5190 16 cm
AURA-7-5191 19 cm
Koser-Konig



AURA-7-5192
13 cm



AURA-7-5193
14.5 cm
Oakland



AURA-7-5194
Maunder Wood
9.5 cm




AURA-7-5195
Spreading 11 cm
47 mm



AURA-7-5196
13 cm



AURA-7-5197 Small
AURA-7-5198 Medium
AURA-7-5199 Large
tor



72 x 23 mm

AURA-7-5200
15 cm



Doy AURA-7-5201
13.5 cm




AURA-7-5202



AURA-7-5203




AURA-7-5204
19 cm




85 x 24 mm

AURA-7-5205
15 cm



88 x 25 mm

AURA-7-5206
16 cm



30 mm

36 mm

Wieder

AURA-7-5207

AURA-7-5208

Lip and Cheek Retractors



AURA-7-5209
hin
Retractor
16.5 cm



16 cm
AURA-7-5210 8 mm
AURA-7-5211 10 mm



AURA-7-5212 43 x 11 mm
AURA-7-5213 57 x 11 mm
AURA-7-5214 83 x 16 mm



35 x 11 mm

AURA-7-521
Obwegeser
Epithel retractor
22 cm



72 x 11 mm

AURA-7-5216
Ramus retractor for
separation of ramus & tissue
22 cm



AURA-7-5217
10.5 cm



30x19 mm

33x25 mm

AURA-7-5218
15 cm



AURA-7-5219
14 cm



AURA-7-5220
esota
15 cm



AURA-7-5221



AURA-7-5222

Fits any simple stem mouth mirror handle



AURA-7-5223

with handle



AURA-7-5224

20 cm

Finger protectors



AURA-7-5225 Small
AURA-7-5226


AURA-7-5227 Large

Langenbeck



AURA-7-5228

11.5 cm



AURA-7-5229

12cm

cheek Retractor



AURA-7-5230

14.5 cm



Sternberg


AURA-7-5231 12 cm
AURA-7-5232 14 cm
AURA-7-5233 15 cm



AURA-7-5234




AURA-7-5235



AURA-7-5236

Implants Retractors


Retractors



Senn-Miller
16 cm

AURA-7-5237

AURA-7-5238



Max
16 cm

AURA-7-5239

AURA-7-5240

AURA-7-5241

AURA-7-5242




Max
16 cm

AURA-7-5243

AURA-7-5244

AURA-7-5245

AURA-7-5246



Kirkland
Gum & Tissue
Retractors

AURA-7-5247

AURA-7-5248

AURA-7-5249

AURA-7-5250




AURA-7-5251
holding forcep
17 cm



e holding forcep
AURA-7-5252 15 cm
AURA-7-5253 17 cm

AURA-7-5254
air




AURA-7-5255
Pair
cm


AURA-7-5256
Pair
nbeck 21 cm


Implantology Instruments

Osteotomes






Convex Tips





Osteotomes with Stop




Implantology Instrumentns

PDL Elevators



AURA-7-5345 Distal, 2.5 mm Micro Edge



AURA-7-5346 Mesial, 2.5 mm Micro Edge



AURA-7-5347 Distal, 4 mm Micro Edge



AURA-7-5348 Mesial 4 mm Micro Edge



AURA-7-5349 2.5 mm Micro Edge, Straight



AURA-7-5350 2.5 mm Micro Edge, Curved



AURA-7-5351 3.4 mm

Straight AURA-7-5352 Set


Of 7 Pieces



AURA-7-5353 Distal, 2.5 mm Micro Edge



AURA-7-5354 Mesial, 2.5 mm Micro Edge



AURA-7-5355 Distal, 4 mm Micro Edge



AURA-7-5356 Mesial 4 mm Micro Edge



AURA-7-5357 2.5 mm Micro Edge, Straight



AURA-7-5358 2.5 mm Micro Edge, Curved



AURA-7-5359 3.4 mm Straight

AURA-7-5360 Set Of 7 Pieces



Instruments Tray

Trepine Drills



AURA-7-5361
Exterior 3 mm



AURA-7-5362
Exterior 4 mm



AURA-7-5363
Exterior 5 mm



AURA-7-5364
Exterior 6 mm



AURA-7-5365
Exterior 7 mm



AURA-7-5366
Exterior 8 mm



AURA-7-5367
Exterior 9 mm



AURA-7-5368
Exterior 10 mm




AURA-7-5369
Drills Stand
Complete with drills



AURA-7-5370
Drills Stand
Complete with drills

Tissue Punches



Keyes



AURA-7-5371



AURA-7-5372



AURA-7-5373



AURA-7-5374



AURA-7-5375



AURA-7-5376



AURA-7-5377



AURA-7-5378



AURA-7-5379



AURA-7-5380



AURA-7-5381



Keyes

Implantology Instrumetns



AURA-7-5382
r Chisel Straight



AURA-7-5383
ight bone collector chisel



AURA-7-5384
r Chisel Curved



AURA-7-5385
ved bone collector chisel



AURA-7-5386
w driver handle
Hexagonal head



AURA-7-5387
zer with Titanium Grinding



AURA-7-5388
Basin
Dia 35 mm



AURA-7-5389
Bone Mill



AURA-7-5390
Titanium

Page 152–155

▪ Impression Trays



All Impression Trays Upper & Lower Type
are available in S, M, L & XL Sizes



AURA-7-5391
"CK"
stainless steel
water-cooled
impression trays
with retention rim.



AURA-7-5392
"CK"
stainless steel
water-cooled impression
trays with retention rim.



AURA-7-5393
"R"
stainless steel
impression trays
with retention rim.




AURA-7-5394
"R" stainless
steel impression trays
with retention rim.



AURA-7-5395, "solid
regular stainless
steel impression
trays without
retention rim.



AURA-7-5396, "solid regular
stainless steel impression
trays without retention rim.



AURA-7-5397
total denture
stainless steel
impression trays with rim.



AURA-7-5398
total denture
stainless steel
impression trays with rim.

Stainless Steel Impression Trays

All Impression Trays Upper & Lower Type are available in S, M, L & XL Sizes



AURA-7-5399 CK" solid regular stainless steel impression trays with retention rim.



AURA-7-5400 CK" solid regular stainless steel impression trays with retention rim.



AURA-7-5401



AURA-7-5402



AURA-7-5403



AURA-7-5404



AURA-7-5405



AURA-7-5406



AURA-7-5407

PACKINGS



AURA-7-5408



AURA-7-5409 ss steel anterior depressed impression trays with retention rim.



AURA-7-5410



AURA-7-5411



AURA-7-5412



AURA-7-5413



AURA-7-5414



AURA-7-5415



AURA-7-5416 ss steel anterior depressed impression trays with retention rim.

All Impression Trays Upper & Lower Type
are available in S, M, L & XL Sizes



AURA-7-5417
ambini - for children



AURA-7-5418
toothed upper jaws



AURA-7-5419
ambini - for children



AURA-7-5420
toothed upper jaws



AURA-7-5421
toothed lower jaws



AURA-7-5422
ous upper jaws



AURA-7-5423
ous lower jaws

Stainless Steel Impression Trays

All Impression Trays Upper & Lower Type are available in S, M, L & XL Sizes



AURA-7-5424 Entire stainless steel impression tray with rim.



AURA-7-5425 Entire stainless steel impression tray with rim.



AURA-7-5426 CK" perforated regular stainless steel impression trays with retention rim.



AURA-7-5427 CK" perforated regular stainless steel impression trays with retention rim.



AURA-7-5428 Stainless steel anterior depressed impression trays with retention rim.



AURA-7-5429



AURA-7-5430



AURA-7-5431



AURA-7-5432



AURA-7-5433



AURA-7-5434



AURA-7-5435 Stainless steel anterior depressed impression trays with retention rim.



AURA-7-5436 Stainless steel impression trays without retention rim for complete denture.



AURA-7-5437 Stainless steel impression trays without retention rim for complete denture.

Page 156–157

- Arriculators



AURA-8-5101



AURA-8-5102



AURA-8-5103



AURA-8-5104



AURA-8-5105



AURA-8-5106

Articulators



AURA-8-5107



AURA-8-5108



AURA-8-5109



AURA-8-5110

Page 158–160

- Dental Surgical Hollowares





AURA-8-5111
ENDODONTICBOX



AURA-8-5112
STERILIZINGBOX
Available in dif. sizes
17x17x3 cm
20x10x4 cm
25x12x6 cm
30x14x6 cm



AURA-8-5113
ENDO TRAY BOX



AURA-8-5114
COTTONDISPENSER



AURA-8-5115
ALCOHOLLAMP

Sterilizing Drums



AURA-8-5116



AURA-8-5117



AURA-8-5118

Dental & Surgical Hollowares

Sterilizing Boxes



AURA-8-5119



AURA-8-5120



AURA-8-5121

Sterilizing Boxes



AURA-8-5122



AURA-8-5123



AURA-8-5124



AURA-8-5125



AURA-8-5126



AURA-8-5127

FORCEP AND THERMOMETER JARS

COTTON DISPENSER



AURA-8-5128
IRRIGATOR(Douchcan)



AURA-8-5129
MEASURINGJUG



AURA-8-5130
KIDNEY TRAY



AURA-8-5131
FORCEP JAR



AURA-8-5132
GALLIPOT WITH & WITHOUT LID



AURA-8-5133
SPITTEN MUG



AURA-8-5134
SPIRIT LAMP



AURA-8-5135
Perforated Tray



AURA-8-5136
Perforated Base



AURA-8-5137
Solid Base



AURA-8-5138
Instruments Tray



AURA-8-5139
Instruments Tray



AURA-8-5140
**Instruments Tray
(Perforated)**



AURA-8-5141
**Instruments Tray Shallow
(Mayo)**



AURA-8-5142
Scalers Tray



AURA-8-5143
Cathetar Tray

Alphabetical Index

A		
Adams	125, 133	
Aderer	133	
Adson	96	
Ainsworth	110	
Albee	77	
Albrecht	109	
Alcohol Lamp	142	
Alginate	93-94	
Allen	109	
Amalgam Carriers	117	
Amalgam Guns	117	
Amalgam Well	117	
Andresen	134	
Angle	133	
Antkoviak	78	
Apical	36, 44	
Arruga	97	
Articulators	140	
Ash	53, 88	
Aspirating Syringes	102	
B		
Band Adaptor	102	
Barry	47, 50	
Baumgarthner	96	
Beale	76	
Beebee	116	
Bein	33, 40	
Bennet	53, 123	
Beyer	121	
Bicuspids	114	
Birdbeak Pliers	132	
Blumenthal	120, 121	
Bohler	120	
Bohm	118	
Boley	77	
Bone Chisels	68	
Boynton	97	
Breech Loading Syringe	102	
Brewer	110	
Buck	61-63, 65	
Buckley	67, 124	
Burnishers	58	
C		
Caroll	77	
Castroviejo	77, 98, 100	
Cathetar Tray	144	
Cement Spatulas	91, 92	
Chompret	44	
Ck1	108	
CK2	108	
Cleveland	122	
College	105, 107	
Columbia	64, 66, 74	
Columbia	79	
Cone-Socket	90	
Cork Screw Tip	104	
Cotton Dispenser	142, 143	
Coupland	37	
Crane-kaplan	74	
Crile-Wood	96	
Crown Removers	115-116	
Cryer	35, 37, 48	
Curette	64, 66	
D		
De La Rosa	129	
Dental Probes	57	
Dental Syringes	103-104	
Derf	95	
Didner	113	
E		
Eliselberg-Mathieu	97	
Elliott	53	
Endo Finger Ruler	102	
Endo Tray Box	142	
Endodontic Box	142	
Ermert	78	
Excavators	87-88	
Explorers	81	
Extracting Forceps	01-31	
F		
Fahnenstock	79	
Filling Instruments	52-56	
Fischer	134	
Flagg	107	
Flohr	38	
Forcepps Jar	143, 144	
Fox	58	
Fraizer	100	
Freer	123	
Friedman	119-120	
Furrer	119	
G		
Gallipot	144	
Gardner	67, 122	
GF1	108	
GF2	108	
Gillies	96	
Goldman Fox	60, 63, 65, 73, 120	
Goslee	134	
Gouges	68	
Gracey Curettes	69	
Gracey Mesial Distal	69	
Gracey	56-69-72	
Gritman	76-78-80	
H		
Halstead-Mosquito	118	
Hegar-Baumgartner	95-96	
Heidbrink	45-48-50	
Heidemann	55	
Hemingway	75	
Henahan	123	
How	133	
Hylin	40, 76	
I		
Impression Trays	136-139	
Intra-Ligamental Syringes	104	
Irrigator	143	
Ivory	110, 113	
Iwanson	77	
J		
Jaquette	74	
Johnson	133	
K		
Kelly	118	
Kerrison	120	
Kidney Tray	143	
Kilian	124	
Kilian-Claus	124	
Klein	07	
Kleinert-Kutz	121	
Kopp	39, 49, 108	
L		
Landgraf	107	
Langenbeck	134	
Le Cron	76	
Lecluse	43	
Lessmann	79	
Lichtenberg	97	
Listion-mini	122	
London Hospital	41	
London-College	107, 109	
Lucas	75	
Luer	120	
Luks	84	

M			W		
Machtou	83		Wakefield		
Mallets	68		Waldron	67	
Mathieu	99, 100		Waldsachs	133	
Matrix Bands	113, 114		Warwick-James	130, 133	
Matrix Retainers	113		Water Syringe	41	
Mayo-Hegar	95		Weingart	101	
Mc Call	73		White	129, 133	
Mc Kellops	133		Williams	52, 60, 83	
Mead	43, 58, 119		Williger	123	
Measuring Jug	143		Winter	46, 50	
Measuring Plate	84		Woodson	54, 55	
Medicon	78				
Meriam	106, 109		Y		
Merrit A	60		Young	133	
Mershon	110, 116		Younger-Good	73	
Micro-Barraquer	99-100				
Micro-Mosquito	118		Z		
Miller	118, 124		Zahle	76	
Mini-Ryder	95		Zaufal-Jansen	121	
Molar	114		Zerfing-Chisel	64, 66	
Morell	116		Zuricher	77	
Mouth Mirror Handles	89, 90				
Mouth Mirrors	89, 90				
MSY	76, 134				
N					
Nabers	83				
Nance	128, 133				
Niro	122				
O					
Obwegeser	124				
Ochsenbein	61, 62, 67				
Olsen-Heger	96				
Oms	83				
Ott	41				
P					
P.K. Thomas	76				
Parts	68				
Parts	134				
Pocket Probes	58-60				
Plaster Knives	80				
Plaster Spatulas	93-94				
Position Gauge	102				
Pott	42				
R					
Rein	81				
Relax	05				
Reynolds	129				
Rhodes	61, 62				
Rhodium	89				
Roach	76				
Root Elevators	32--51				
Root Splinter	31				
Root-Canal Excavators	86				
Root-Canal Explorers	86				
Root-Canal Plungers	84-85				
Rubber Dam Clamps	111-112				
Rubber Dam Frame	110				
Ruskin-Mini	121				
S					
Scalers	58				
Schell	119				
Schilder	85				
Schluger	61, 62				
Schwarz	134				
Seldin	34, 48, 123				
Spirit Lamp	144				
Spitten Mug	144				
Spreader	82-83				
Sterilizing Boxes	142, 143				
Sterilizing Drums	142				
Sugarman	61, 162				
T					
Taylor	74				
Toennis	97				
Tofflemire	114				
Tofflemire-Junior	113				
Tofflemire-Senior	113				
Tofflemire-Universal	113				
Towner Jaquette	74				
Towner	74				
Trays	144				
Treymann	116				
Tweed	130				
U					
Univarsal "HM"	116				
Univarsal	128, 133				
University of Michigan	60				
USA Model	133				
V					
Vehe	76				
Vernier	77				
Vertical Condenser	83				